



LymphoTrack[®] Assays


Technical Training

MiSeq[®] and Ion S5/PGM[™] Platforms

LymphoTrack[®] - Workflow



LymphoTrack[®] - MiSeq[®] Target Loci



LymphoTrack[®] - S5/PGM[™] Target Loci



IGH
FR1

IGH
FR2

IGH
FR3

IGK

TRG



LymphoTrack® - Available Kits

RUO Products	MiSeq®		S5/PGM™ 12 indices (60 rxn)
	Kit A 8 indices (40 rxn)	Panel 24 indices (120 rxn)	
LymphoTrack® IGHV Leader Somatic Hypermutation Assay	✓	✓	
LymphoTrack® IGH FR1 Assay	✓	✓	✓
LymphoTrack® IGH FR2 Assay	✓	✓	✓
LymphoTrack® IGH FR3 Assay	✓	✓	✓
LymphoTrack® IGH FR1/2/3 Assay	✓	✓	✓
LymphoTrack® IGK Assay	✓	✓	✓
LymphoTrack® TRG Assay	✓	✓	✓
LymphoTrack® TRB Assay	✓	✓	

48 indices available for IGH FR1 MiSeq®

LymphoTrack[®] Assays

Experiment Planning




LymphoTrack[®] - Experiment Planning

Key Factors

- How many samples?
- How many targets?
- Which targets?

These will determine:

- # of indices needed (8 / 12 / 24)
- MiSeq[®] Kit/Cycles, Ion[™] Chip
- Cost per target per sample



LymphoTrack[®] - MiSeq[®] Compatibility

LymphoTrack [®] Assay	Avg. Target Size incl. Target, Index & Adaptor (bp)	MiSeq [®] Reagent		
		Kit v2 (300 cycle)	Kit v2 (500 cycle)	Kit v3 (600 cycle)
IGHV Leader SHM	660	X	X	○
IGH FR1	450	X	○	○
IGH FR2	390	X	○	○
IGH FR3	260	○	○	○
IGK	410	X	○	○
TRG	300	○	○	○
TRB	400	X	○	○

MiSeq [®] Reagent Kit	Read Length	MCS Version	MiSeq [®] Run Time	Reads	%≥Q30
v2	2 x 151 bp	v2.6	~ 24 hours	up to 15 Million	> 80%
	2 x 251 bp	v2.6	~ 39 hours		> 75%
v3	2 x 301 bp	v2.6	~ 56 hours	up to 25 million	> 70%

LymphoTrack® - S5/PGM™ Compatibility

LymphoTrack® Assay	Avg. Target Size incl. Target, Index & Adaptor (bp)	Chef/OT2™ Reagent for S5/PGM™
IGH FR1	430	<p><u>For Ion S5™ use:</u> Ion 510, Ion 520, or Ion 530 – for Ion Chef™ or Ion 520 Chip or Ion 530 Chip - for OneTouch 2™</p> <p><u>For Ion PGM™ use:</u> PGM Hi-Q View</p>
IGH FR2	370	
IGH FR3	240	
IGK	390	
TRG	280	

Platform	Template Prep.	Sequencing Kit	Chip	Reads
Ion S5™	Ion Chef™	Ion 510 – Chef Ion 520 – Chef Ion 530 – Chef	Ion 520 Chip	4 – 6 million
	OneTouch 2™	Ion 520 – OT2 Ion 530 – OT2	Ion 530 Chip	15 – 20 million
Ion PGM™	OneTouch 2™	PGM Hi-Q View – OT2	Ion 316 Chip kit v2 BC	2 – 3 million
			Ion 318 Chip Kit v2 BC	4 - 5.5 million

LymphoTrack® - Multiplexing

Single Assay Multiplexing

- Up to 12 (S5/PGM™) or 24 (MiSeq®) samples and controls can be multiplexed into one run
- Up to 48 indices available for IGH FR1 MiSeq®

Multi-Target Multiplexing


- The same sample can be tested across multiple LymphoTrack® assays in one run
- Each sample receives a unique index
- # of sequencing cycles must be sufficient for the largest amplicon in the multiplex

	1	2	3	4
A	Index 01 Sample A	Index 09 Sample I	Index 01 Sample A	Index 09 Sample I
B	Index 02 Sample B	Index 10 Sample J	Index 02 Sample B	Index 10 Sample J
C	Index 03 Sample C	Index 11 Sample K	Index 03 Sample C	Index 11 Sample K
D	Index 04 Sample D	Index 12 Sample L	Index 04 Sample D	Index 12 Sample L
E	Index 05 Sample E	Index 13 Sample M	Index 05 Sample E	Index 13 Sample M
F	Index 06 Sample F	Index 14 PosCtrl	Index 06 Sample F	Index 14 PosCtrl
G	Index 07 Sample G	Index 15 NegCtrl	Index 07 Sample G	Index 15 NegCtrl
H	Index 08 Sample H	Index 15 NTC	Index 08 Sample H	Index 15 NTC
	IGH FR1		TRG	

LymphoTrack® - Example Calculation

A MiSeq® v3 Flow Cell outputs >20 million reads

$$\begin{array}{r} 20\,000\,000 \text{ Reads} \\ \div 24 \text{ Samples} \\ \hline \div 7 \text{ Targets} \\ \hline \sim 120\,000 \text{ Reads per} \\ \text{Target, Per Sample} \end{array}$$



LymphoTrack[®] Assays

Workflow Overview




LymphoTrack[®] - Workflow

All LymphoTrack[®] Assays follow the same **7 step procedure**


Multiplexed assays can be run side-by-side to **reduce laboratory time**

This simplified process allows even new NGS users to succeed with ease




LymphoTrack[®] - MiSeq[®] Time

<u>Step</u>	<u>Hands-On</u>	<u>Total Time</u>
Prepare DNA	30 min	30 min
Run PCR	0 min	150 min
Purify Amplicons	30 min	45 min
Quantify	60 min	180 min
Dilute & Pool	30 min	30 min
Load Instrument	30 min	30 min
Sequencing	0 min	~48 hours
Total Time	~3 hours	~56 hours



LymphoTrack[®] - S5/PGM[™] Time

<u>Step</u>	<u>Hands-On</u>	<u>Total Time</u>
Prepare DNA	30 min	30 min
Run PCR	0 min	150 min
Purify Amplicons	30 min	45 min
Quantify	15 min	45 min
Dilute & Pool	30 min	30 min
Load Instrument	OT2: 120 min Chef: 30 min	OT2: 7 hr Chef: 12 hr
Sequencing	30 min	~4 hours
Total Time	~3-5 hours	~20-25 hours



Ion S5/PGM™ - Technical Differences

LymphoTrack[®] Procedures are consistent across both platforms

Ion Chef™ helps automate template preparation

- Minimizes hands-on time
- Improves consistency of results



Ion S5™ is compatible with Ion OT2™, Chef™ & LymphoTrack[®] Assays

- Ability to generate longer amplicons
- Increased number of reads



Ion Chef™ & Ion S5™ provide a simple to use, automated NGS solution

- Cartridge-based reagents
- Consumable tracking with the automated RFID

Ion S5/PGM™ - Run Time Comparisons

Step	Total Processing Time (Hours)		
	OT2/PGM	OT2/S5	Chef/S5
PCR	~2.5	~2.5	~2.5
AMPure XP Purification	~0.5	~0.5	~0.5
Quantification & Pooling	~1.0	~1.0	~1.0
Template Prep	~6.5*	~6.5*	~12.0*
Sequencing*	~4.0 – 8.0	~4.0	~4.0
Time-to-Result	~14.5 – 18.5	~14.5	~20.0

Sequencing time will vary depending on the chip in use
 Analysis time will vary depending on instrument and chip type

*Hands-on time is significantly reduced for Ion Chef

Step 1 - Prepare DNA Samples

DNA Input

- Use an appropriate resuspension/ elution buffer such as 1/10 TE buffer (1 mM Tris-HCl, 0.1 mM EDTA, pH 8.0)
- Quantify DNA with a method specific for double-stranded DNA (dsDNA) e.g. Qubit.
- Minimum input quantity is 50 ng of high-quality DNA in 5 μ L volume (10 ng/ μ L)
- Invivoscribe Specimen Control Size Ladder can be used to confirm quality of DNA.



Step 2 - Run PCR

Simplified Setup

- Premade master mix: **just add enzyme and sample**
- Always include positive, negative, and non-template controls!

Reagent	Volume
Master Mix	45 µL
EagleTaq™ DNA polymerase	0.2 µL
Sample or Control DNA	5 µL
Total Volume	50.2 µL

Step	Temperature	Time	Cycles
1	95 °C	07:00	1
2	95 °C	00:45	29x*
3	60 °C	00:45	
4	72 °C	01:30	
5	72 °C	10:00	1
6	15 °C	Hold	1

*32x for IGHV Leader



Step 3 - Purify Amplicons

SPRI Bead Purification

- Removes excess primers, salts, and enzymes

Purification Procedure

- Add beads to each sample and incubate at room temperature
- Place on magnetic stand to separate beads from the supernatant
- Wash twice with 200 μ L of fresh 70/80% Ethanol
- Air dry at room temperature
- Elute in 10mM Tris pH 8.0



Step 4 - Quantification

PCR Product Quantification (MiSeq[®])

- Quantify PCR products using KAPA Library Quantification Kit for Illumina systems
- Compare Ct values of positive and negative controls to the NTC (ΔCt must be $\geq 4.0^*$)

PCR Product Quantification (S5/PGM[™])

- Quantify PCR products using either 2100 BioAnalyzer or LabChip Gx systems

Library Pooling & Quantification

- Dilute samples to 4 nM in at least 10 μ L and then pool them together to be equimolar
- Quantify library using the same quantification platform
- MiSeq[®]: Denature, and then dilute pool to 12-20 pM
- S5/PGM[™]: Dilute pool to 20 pM


* ΔCt must be ≥ 3.0 for TRB MiSeq Assay

Optional - Target Multiplexing

Generate sample specific library pools & quantify

Use target size and pool concentrations to create an equimolar final loading pool

Best practice to quantify final loading pool for accurate loading



Step 5 - Load and Run MiSeq®

Library Denaturation & Dilution

- The library pool is denatured into single strands using sodium hydroxide
- After incubation, the library is diluted to 40 pM

Loading Concentration

- One final dilution prepares the library for loading.
- The final loading concentration is dependent on the MiSeq® chemistry being used and the targets being sequenced
- Libraries should be loaded within an hour of denaturing & diluting

Reagent	Volume
4 nM Library	10 µL
0.2 N NaOH	10 µL
Incubate for 5 minutes at RT	
Cold HT1 Buffer	980 µL
Total (40 pM)	1000 µL

Reagent	Loading Concentration		
	12 pM	15 pM	20 pM
40 pM Library	300 µL	375 µL	500 µL
HT1 Buffer	700 µL	625 µL	500 µL
TOTAL	1000 µL	1000 µL	1000 µL

Step 5 - Load and Run MiSeq®

Loading

- Load 600 µL of the Final Prepared Library onto a MiSeq[®] Reagent Cartridge

Sample Sheet


- Use Illumina Experiment Manager™ OR
- Use the provided SampleSheet.csv file and modify the template for your sample needs

Starting Your Run

- Direct the MiSeq[®] to the saved sample sheet
- Start the instrument using the user interface



Step 5 - Load and Run S5/PGM™




Sequencing by Synthesis


A quick look at what is happening
inside the sequencer



Sequencing by Synthesis




Sequencing by Synthesis



MiSeq® Base Calling

Fluorescence Detection

- Nucleotides with fluorophores are used to extend from the sequencing primer
- The added base is excited with light and the resulting fluorescence is imaged by a camera
- The color that is emitted allows the instrument to determine which base was added
- The dye and terminator are removed and the process is repeated for the next cycle/base



<https://binf.snipcademy.com/lessons/ngs-techniques/illumina-solexa>

Reads

- When a single strand of DNA is sequenced, this is considered one '**read**'
- During PCR, many copies of each starting DNA strand were generated
- Copies of a normal rearrangement can be sequenced, producing **multiple reads with the same sequence**
- When a very large number of reads have the same sequence, it likely means there were more than one of the original strand
- This is the **basis for clonality testing by NGS**: an **overabundance in reads** of one sequence when compared to a background of relatively low prevalence (polyclonal) reads

<https://binf.snipcademy.com/lessons/ngs-techniques/illumina-solexa>

LymphoTrack[®] Assays

Workflow Overview (continued)




Step 6 - Included Bioinformatics

Data Analysis

- Transfer the appropriate files to a local folder (not on the Instrument)
- Launch LymphoTrack[®] Software
- Select targets, decimal format, and location of .fastq files
- Select 'Launch Program'

LymphoTrack[®] MiSeq Data Analysis version 2.4.3 Copyright 2018 Invivoscribe Inc.
FOR RESEARCH USE ONLY; not for use in diagnostic procedures.

 **LymphoTrack[®]**
Software

Choose analysis mode(s):

Leader TRG TRB IGH FR.1 IGH FR.2 IGH FR.3 IGH

Choose decimal mark:

"." decimal point, e.g., "1,000.00"

"," decimal comma, e.g., "1.000,00"

Choose an Input Folder for analysis:

Step 6 - Included Bioinformatics

PDF reports are automatically generated per sample and target

Data is presented numerically to allow for **objective interpretation**

Graphs are provided to allow for visualization of results


Total Read Count: 673441

IndexQ30: 89.45

Caution: Do not edit fields and save.

Top 10 Merged Read Summary

Rank	Sequence	Length	Merge count	V-gene	J-gene	% total reads	Cumulative %	Mutation rate to partial V-gene (%)	In-frame (Y/N)	No Stop codon (Y/N)	V-coverage	CDR3 Seq
1	CATCTGGATACAC ⁺	295	49206	IGHV1-46_03	IGHJ4_02	7.31	7.31	0.00	Y	Y	100.00	not found
2	GCCTCTGGATTCA ⁺	300	181	IGHV3-35_01	IGHJ6_02	0.03	7.33	3.96	N	N	99.12	not found
3	ACCTCTGCAATCA ⁺	159	172	IGHV3-21_04	IGHJ6_02	0.03	7.36	1.32	n/a	N	22.03	not found
4	ACCTCTGGATTTC ⁺	272	171	IGHV3-74_02	IGHJ4_02	0.03	7.38	6.67	Y	Y	88.44	not found
5	GCGTCTGGAATCA ⁺	166	160	IGHV3-30_02	none	0.02	7.41	5.73	n/a	N	72.25	not found



Step 7 - Analyze

MiSeq[®] Validity Specs

- MiSeq[®] instrument run validity specifications are dependent on the MiSeq[®] kit chemistry and the number of sequencing cycles

S5/PGM[™] Run Validity Specs

- Loading > 50%
- Enrichment > 50%
- Clonal > 50%

Assay Specifications

- Positive Control top % Reads $\geq 2.5\%$
- Negative Control top % Reads $< 1.0\%$
- If all validity criteria are met, move forward to interpretation & reporting

Step 7 – Analyze

Top 10 Merged Read Summary

- The Read Summary tab only shows the top 10 reads after merging with the top 500 reads that differ in 1 or 2 nucleotides

Note: This analysis is also possible with *IGH FR1*

LymphoTrack Report for assay LEADER

Sample name: Leader_positive_S23_L001_001_combined

Total Read Count: 474947

IndexQ30: 87.88

Caution: Do not edit fields and save.

Top 10 Merged Read Summary

Rank	Sequence	Length	Merge count	V-gene	J-gene	% total reads	Cumulative %	Mutation rate to partial V-gene (%)	In-frame (Y/N)	No Stop codon (Y/N)	V-coverage	CDR3 Seq
1	TTCTCGTGGTGGC	455	50248	IGHV4-59_08	IGHJ4_02	10.58	10.58	11.26	Y	Y	98.63	GCGAGACGGAGC
2	CTGCTACTGACTG	319	192	IGHV2-70_10	IGHJ4_02	0.04	10.62	4.32	n/a	N	35.55	not found
3	CTGCTGCTGACCA	466	175	IGHV2-5_01	IGHJ5_01	0.04	10.66	6.62	Y	Y	100.00	GCACACAGACCGC
4	CTGCTGCTGACCA	457	162	IGHV2-5_05	IGHJ6_02	0.03	10.69	2.99	Y	Y	99.67	GCACACAGATACT
5	CTGCTGCTGACCA	474	154	IGHV2-5_05	IGHJ4_02	0.03	10.72	3.99	Y	Y	99.67	GCACACAGATACT
6	CTGCTGCTGACCA	454	150	IGHV2-5_10	IGHJ5_02	0.03	10.76	11.78	Y	Y	98.99	GCATATGGTGTA
7	CTGCTGCTGACCA	469	139	IGHV2-5_01	IGHJ4_02	0.03	10.78	1.32	Y	Y	97.68	GCACTCGCGACAC
8	CTCGCCCTCCTCC	466	139	IGHV5-51_01	IGHJ4_02	0.03	10.81	7.09	Y	Y	99.32	GCGAGATACTATT
9	CTGCTACTGACTG	490	137	IGHV2-70_10	IGHJ3_02	0.03	10.84	0.66	Y	Y	99.34	GCACGGATTCTC
10	CTGCTGCTGACCA	478	135	IGHV2-5_10	IGHJ6_02	0.03	10.87	3.70	Y	Y	98.99	GCATACACTTGT

Step 7 – Analyze

LymphoTrack Report for assay LEADER

Sample name: Leader_positive_S23_L001_001_combined

Total Read Count: 474947

IndexQ30: 87.88

Caution: Do not edit fields and save.

Top 10 Merged Read Summary

Rank	Sequence	Length	Merge count	V-gene	J-gene	% total reads	Cumulative %	Mutation rate to partial V-gene (%)	In-frame (Y/N)	No Stop codon (Y/N)	V-coverage	CDR3 Seq
1	TTCTCGTGGTGGC+	455	50248	IGHV4-59_08	IGHJ4_02	10.58	10.58	11.26	Y	Y	98.63	GCGAGACGGAGC+
2	CTGCTACTGACTG+	319	192	IGHV2-70_10	IGHJ4_02	0.04	10.62	4.32	n/a	N	35.55	not found
3	CTGCTGCTGACCA+	466	175	IGHV2-5_01	IGHJ5_01	0.04	10.66	6.62	Y	Y	100.00	GCACACAGACCGC+
4	CTGCTGCTGACCA+	457	162	IGHV2-5_05	IGHJ6_02	0.03	10.69	2.99	Y	Y	99.67	GCACACAGATACT+
5	CTGCTGCTGACCA+	474	154	IGHV2-5_05	IGHJ4_02	0.03	10.72	3.99	Y	Y	99.67	GCACACAGATACT+
6	CTGCTGCTGACCA+	454	150	IGHV2-5_10	IGHJ5_02	0.03	10.76	11.78	Y	Y	98.99	GCATATGGTGTA+
7	CTGCTGCTGACCA+	469	139	IGHV2-5_01	IGHJ4_02	0.03	10.78	1.32	Y	Y	97.68	GCACTCGCGACAC+
8	CTCGCCCTCCTCC+	466	139	IGHV5-51_01	IGHJ4_02	0.03	10.81	7.09	Y	Y	99.32	GCGAGATACTATT+
9	CTGCTACTGACTG+	490	137	IGHV2-70_10	IGHJ3_02	0.03	10.84	0.66	Y	Y	99.34	GCACGGATTCCCTG+
10	CTGCTGCTGACCA+	478	135	IGHV2-5_10	IGHJ6_02	0.03	10.87	3.70	Y	Y	98.99	GCATACACTTGT+

Expected Values IGHV SHM

- IGH Positive Control top % reads $\geq 2.5\%$
- NGS Negative Control top % reads $< 1.0\%$
- IGH SHM Positive Control (4-088-0008, can be purchased separately) top % reads $\geq 2.5\%$
- IGH SHM Positive Control mutation rate $\geq 2.0\%$
- MiSeq Run Validity Q30 $> 70\%$ for v3 (2x301)

Optional - Raw Data

Merged Top 10 File


- 10 most prevalent merged sequences
- Top 500 reads are scanned for having up to 2 mismatches to the top 10 sequences to account for PCR & Sequencing error

Read Summary File

- Top 200 unmerged reads
- Aligned to identified V and J genes

Unique Reads File

- All unique reads identified in the fastq file. Can be useful for troubleshooting



Take Home Message

All LymphoTrack[®] assays follow **one simple workflow**

NGS can seem overwhelming – **LymphoTrack[®] simplifies the process** by reducing steps, labor costs, and analysis time

Multiplexing targets into one run is a great way to **reduce cost per target**

Included bioinformatics package **streamlines analysis** and provides flexible reporting methods